



# CRD1 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA862988LD01DOA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9M591
<b>Immunogen</b>	Recombinant Arabidopsis thaliana Magnesium-protoporphyrin IX monomethyl ester [oxidative] cyclase, chloroplastic protein (37-409AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Arabidopsis thaliana
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Magnesium-protoporphyrin IX monomethyl ester [oxidative] cyclase, chloroplastic (Mg-protoporphyrin IX monomethyl ester oxidative cyclase) (EC 1.14.13.81) (Copper response defect 1 protein) (Dicarboxylate diiron protein) (AtZIP) (MPE-cyclase), CRD1, ACSF AT103 CHL27 ZIP
<b>Immunogen Species</b>	Arabidopsis thaliana
<b>Research Area</b>	Others
<b>Target Names</b>	CRD1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.