



VDAC1 Antibody, Biotin conjugated

Product Code	CSB-PA865954LD01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9SRH5
Immunogen	Recombinant Arabidopsis thaliana Mitochondrial outer membrane protein porin 1 protein (2-276AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial outer membrane protein porin 1 (Voltage-dependent anion-selective channel protein 1) (AtVDAC1) (VDAC-1), VDAC1
Immunogen Species	Arabidopsis thaliana
Research Area	Tags & Cell Markers
Target Names	VDAC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.