



PEX14 Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA867033LB01DOA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9FXT6 |
| Immunogen | Recombinant Arabidopsis thaliana Peroxisomal membrane protein PEX14 protein (288-507AA) |
| Raised In | Rabbit |
| Species Reactivity | Arabidopsis thaliana |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Peroxisomal membrane protein PEX14 (Peroxin-14) (AtPEX14) (Peroxisome biogenesis protein 14) (Pex14p) (Protein PEROXISOME DEFECTIVE 2), PEX14, PED2 |
| Immunogen Species | Arabidopsis thaliana |
| Research Area | Tags & Cell Markers |
| Target Names | PEX14 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |