



PEX14 Antibody, Biotin conjugated

Product Code	CSB-PA867033LD01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9FXT6
Immunogen	Recombinant Arabidopsis thaliana Peroxisomal membrane protein PEX14 protein (288-507AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisomal membrane protein PEX14 (Peroxin-14) (AtPEX14) (Peroxisome biogenesis protein 14) (Pex14p) (Protein PEROXISOME DEFECTIVE 2), PEX14, PED2
Immunogen Species	Arabidopsis thaliana
Research Area	Tags & Cell Markers
Target Names	PEX14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.