



C5AR2 Antibody, Biotin conjugated

Product Code	CSB-PA868390PD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P296
Immunogen	Peptide sequence from Human C5a anaphylatoxin chemotactic receptor 2 protein (2-19AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	C5a anaphylatoxin chemotactic receptor 2 (Complement component 5a receptor 2) (G-protein coupled receptor 77), C5AR2, C5L2 GPR77
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	C5AR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.