



C5AR2 Antibody, Biotin conjugated

Product Code	CSB-PA868390PD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P296
Immunogen	Peptide sequence from Human C5a anaphylatoxin chemotactic receptor 2 protein (2-19AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	C5a anaphylatoxin chemotactic receptor 2 (Complement component 5a receptor 2) (G-protein coupled receptor 77), C5AR2, C5L2 GPR77
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	C5AR2