




HMGXB4 Antibody

Product Code	CSB-PA871514PA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGU5
Immunogen	Peptide sequence from Human HMG domain-containing protein 4 protein (489-509AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	HMG domain-containing protein 4 (HMG box-containing protein 4) (High mobility group protein 2-like 1) (Protein HMGBCG), HMGXB4, HMG2L1 HMGBCG
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HMGXB4

Image



Western Blot

Positive WB detected in: Jurkat whole cell lysate, K562 whole cell lysate

All lanes: HMGXB4 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 66 kDa

Observed band size: 66 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.