



HMGXB4 Antibody, FITC conjugated

Product Code	CSB-PA871514PC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGU5
Immunogen	Peptide sequence from Human HMG domain-containing protein 4 protein (489-509AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	HMG domain-containing protein 4 (HMG box-containing protein 4) (High mobility group protein 2-like 1) (Protein HMGBCG), HMGXB4, HMG2L1 HMGBCG
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HMGXB4