



HMGXB4 Antibody, Biotin conjugated

Product Code	CSB-PA871514PD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGU5
Immunogen	Peptide sequence from Human HMG domain-containing protein 4 protein (489-509AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	HMG domain-containing protein 4 (HMG box-containing protein 4) (High mobility group protein 2-like 1) (Protein HMGBCG), HMGXB4, HMG2L1 HMGBCG
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HMGXB4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.