



ZRANB1 Antibody, HRP conjugated

Product Code	CSB-PA883378LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGI0
Immunogen	Recombinant Human Ubiquitin thioesterase ZRANB1 protein (111-397AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin thioesterase ZRANB1 (EC 3.4.19.12) (TRAF-binding domain-containing protein) (hTrabid) (Zinc finger Ran-binding domain-containing protein 1), ZRANB1, TRABID
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZRANB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.