



MED37A Antibody, Biotin conjugated

Product Code	CSB-PA888748LD01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9LKR3
Immunogen	Recombinant Arabidopsis thaliana Mediator of RNA polymerase II transcription subunit 37a protein (28-669AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mediator of RNA polymerase II transcription subunit 37a (Heat shock 70 kDa protein 11) (Heat shock protein 70-11) (AtHsp70-11) (Luminal-binding protein 1) (AtBP1) (BiP1), MED37A, BIP1 HSP70-11 MED37_5
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	MED37A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.