



AXIN2 Antibody, HRP conjugated

Product Code	CSB-PA897472OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2T1
Immunogen	Peptide sequence from Human Axin-2 protein (56-73AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Axin-2 (Axin-like protein) (Axil) (Axis inhibition protein 2) (Conductin), AXIN2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	AXIN2