



AXIN2 Antibody, Biotin conjugated

Product Code	CSB-PA897472OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2T1
Immunogen	Peptide sequence from Human Axin-2 protein (56-73AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Axin-2 (Axin-like protein) (Axil) (Axis inhibition protein 2) (Conductin), AXIN2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	AXIN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.