



SYCE3 Antibody, FITC conjugated

Product Code	CSB-PA003820LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1L190
Immunogen	Recombinant Human Synaptonemal complex central element protein 3 protein (1-88AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptonemal complex central element protein 3 (Testis highly expressed gene 2 protein) (THEG-2), SYCE3, C22orf41 THEG2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SYCE3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.