



ENG Antibody, FITC conjugated

Product Code	CSB-PA007667LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17813
Immunogen	Recombinant Human Endoglin protein (32-259AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoglin (CD antigen CD105), ENG, END
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ENG