



FLCN Antibody, Biotin conjugated

Product Code	CSB-PA008711HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFG4
Immunogen	Recombinant Human Folliculin protein (291-579AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Folliculin (BHD skin lesion fibrofolliculoma protein) (Birt-Hogg-Dube syndrome protein), FLCN, BHD
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	FLCN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.