



NXN Antibody, Biotin conjugated

Product Code	CSB-PA016224HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6DKJ4
Immunogen	Recombinant Human Nucleoredoxin protein (155-435AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoredoxin (EC 1.8.1.8), NXN, NRX
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NXN