



RPL24 Antibody, FITC conjugated

Product Code	CSB-PA020202LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P83731
Immunogen	Recombinant Human 60S ribosomal protein L24 protein (1-157AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L24 (60S ribosomal protein L30) (Large ribosomal subunit protein eL24), RPL24
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPL24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.