



RPP38 Antibody, FITC conjugated

Product Code	CSB-PA020350LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78345
Immunogen	Recombinant Human Ribonuclease P protein subunit p38 protein (210-283AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease P protein subunit p38 (RNaseP protein p38) (EC 3.1.26.5), RPP38
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPP38
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.