



RPP38 Antibody, Biotin conjugated

Product Code	CSB-PA020350LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78345
Immunogen	Recombinant Human Ribonuclease P protein subunit p38 protein (210-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease P protein subunit p38 (RNaseP protein p38) (EC 3.1.26.5), RPP38
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPP38