




UBC Antibody

Product Code	CSB-PA03742A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0CG48
Immunogen	Peptide sequence from Human Polyubiquitin-C protein (19-36AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polyubiquitin-C [Cleaved into: Ubiquitin], UBC
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	UBC

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate, 293T whole cell lysate, HeLa whole cell lysate, MCF-7 whole cell lysate, K562 whole cell lysate, LO2 whole cell lysate, Jurkat whole cell lysate

All lanes: UBC antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 78 kDa

Observed band size: 78 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.