



UBC Antibody, FITC conjugated

Product Code	CSB-PA03742C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0CG48
Immunogen	Peptide sequence from Human Polyubiquitin-C protein (19-36AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polyubiquitin-C [Cleaved into: Ubiquitin], UBC
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	UBC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.