



NT5E Antibody, FITC conjugated

Product Code	CSB-PA10529C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21589
Immunogen	Recombinant Human 5'-nucleotidase protein (32-264AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	5'-nucleotidase (5'-NT) (EC 3.1.3.5) (Ecto-5'-nucleotidase) (CD antigen CD73), NT5E, NT5 NTE
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NT5E
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.