



bepA Antibody, HRP conjugated

Product Code	CSB-PA355426LB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P66948
Immunogen	Recombinant Escherichia coli Beta-barrel assembly-enhancing protease protein (214-487AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-barrel assembly-enhancing protease (EC 3.4.-.-), bepA, yfgC
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	bepA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.