



STXBP5 Antibody, HRP conjugated

Product Code	CSB-PA722568LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T5C0
Immunogen	Recombinant Human Syntaxin-binding protein 5 protein (515-610AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntaxin-binding protein 5 (Lethal(2) giant larvae protein homolog 3) (Tomosyn-1), STXBP5, LLGL3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	STXBP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.