



RASL11B Antibody, HRP conjugated

Product Code	CSB-PA859410LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BPW5
Immunogen	Recombinant Human Ras-like protein family member 11B protein (1-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-like protein family member 11B, RASL11B
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	RASL11B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.