



HIST1H2BL Antibody, Biotin conjugated

Product Code	CSB-PA860357ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99880
Immunogen	Recombinant Human Histone H2B type 1-L protein (2-126AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone H2B type 1-L (Histone H2B.c) (H2B/c), HIST1H2BL, H2BFC
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HIST1H2BL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.