



DUSP15 Antibody, FITC conjugated

Product Code	CSB-PA861994LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1R2
Immunogen	Recombinant Human Dual specificity protein phosphatase 15 protein (1-295AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase 15 (EC 3.1.3.16) (EC 3.1.3.48) (VH1-related member Y) (Vaccinia virus VH1-related dual-specific protein phosphatase Y), DUSP15, C20orf57 VHY
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DUSP15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.