



DUSP15 Antibody, Biotin conjugated

Product Code	CSB-PA861994LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1R2
Immunogen	Recombinant Human Dual specificity protein phosphatase 15 protein (1-295AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase 15 (EC 3.1.3.16) (EC 3.1.3.48) (VH1-related member Y) (Vaccinia virus VH1-related dual-specific protein phosphatase Y), DUSP15, C20orf57 VHY
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DUSP15