



VP40 Antibody, HRP conjugated

Product Code	CSB-PA762349HB01ZAT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q77DJ6
Immunogen	Recombinant Zaire ebolavirus Matrix protein VP40 protein (1-326AA)
Raised In	Rabbit
Species Reactivity	Zaire ebolavirus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrix protein VP40 (Membrane-associated protein VP40), VP40
Immunogen Species	Zaire ebolavirus
Research Area	Others
Target Names	VP40