



FABP4 Recombinant Monoclonal Antibody

Product Code	CSB-RA007945A0HU
Abbreviation	Fatty acid-binding protein, adipocyte
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15090
Immunogen	A synthesized peptide derived from human FABP4
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Lipid transport protein in adipocytes. Binds both long chain fatty acids and retinoic acid. Delivers long-chain fatty acids and retinoic acid to their cognate receptors in the nucleus.
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity-chromatography
Isotype	Rabbit IgG
Clonality	Monoclonal
Alias	Fatty acid-binding protein, adipocyte, Adipocyte lipid-binding protein, ALBP, Adipocyte-type fatty acid-binding protein, A-FABP, AFABP, Fatty acid-binding protein 4, FABP4
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FABP4
Clone No.	1H4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.