



FABP4 Antibody

Product Code	CSB-RA007945A0HU
Abbreviation	Fatty acid-binding protein, adipocyte
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15090
Immunogen	A synthesized peptide derived from human FABP4
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Lipid transport protein in adipocytes. Binds both long chain fatty acids and retinoic acid. Delivers long-chain fatty acids and retinoic acid to their cognate receptors in the nucleus.
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity-chromatography
Isotype	Rabbit IgG
Clonality	Monoclonal
Alias	Fatty acid-binding protein, adipocyte, Adipocyte lipid-binding protein, ALBP, Adipocyte-type fatty acid-binding protein, A-FABP, AFABP, Fatty acid-binding protein 4, FABP4
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Gene Names	FABP4
Accession NO.	1H4