



# Goat Anti-Rabbit IgG(H+L) Antibody; Cy3-conjugated

<b>Product Code</b>	CSB-PA956124GA01Rb
<b>Storage</b>	Store at -20°C. Avoid repeated freeze / thaw cycles.
<b>Raised In</b>	Goat
<b>Species Reactivity</b>	Rabbit
<b>Tested Applications</b>	IF, Flow cyt; Recommended dilution: IF:1:20-1:100
<b>Form</b>	Liquid
<b>Storage Buffer</b>	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6
<b>Purification Method</b>	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Product Type</b>	Secondary Antibody
<b>Immunogen Species</b>	Oryctolagus cuniculus (Rabbit)
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.