



Goat Anti-Rabbit IgG(H+L) Antibody; Cy3-conjugated

Product Code	CSB-PA956124GA01Rb
Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
Raised In	Goat
Species Reactivity	Rabbit
Tested Applications	IF,Flow cyt; Recommended dilution: IF:1:20-1:100
Form	Liquid
Storage Buffer	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6
Purification Method	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Product Type	Secondary Antibody
Immunogen Species	Oryctolagus cuniculus (Rabbit)