




Recombinant Mouse Kit ligand (Kitlg), partial (Active)

Product Code	CSB-AP004071MO
Uniprot No.	P20826
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 μ m filtered 1xPBS, pH 7.4
Product Type	Growth Factor
Immunogen Species	Mus musculus (Mouse)
Biological Activity	The ED50 as determined by the dose-dependent stimulation of TF-1 cells is 4-12 ng/ml.
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	KEICGNPVTDNVKDITKLVANLPNDYMITLNYVAGMDVLP SHCWLRDMVIQLSL SLTTLLDKFSNISEGLSNYSIIDKLGKIVDDL VLCMEENAPKNIKESPKRPETRSF TPEFFSIFNRSIDAFKDFMVASDTSDCVLSSTLGP EKDSRVSVTKPFMLPPVA
Research Area	Immunology
Source	E.coli
Target Names	Kitlg
Expression Region	26-189aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	18.4 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Endotoxin	Less than 0.01 EU/ μ g as determined by LAL method.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.