




Recombinant Human Interleukin-4 receptor subunit alpha (IL4R), partial (Active)

Product Code	CSB-AP004271HU
Abbreviation	Recombinant Human IL4R protein, partial (Active)
Uniprot No.	P24394
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 μ m filtered 1xPBS, pH 7.4
Product Type	Interleukin
Immunogen Species	Homo sapiens (Human)
Biological Activity	The ED50 as determined by its ability to inhibit IL-4-dependent proliferation of TF?1 human erythroleukemic cells is 5-20 ng/ml.
Purity	\geq 95% as determined by SDS-PAGE.
Sequence	MKVLQEPTCVSDYMSISTCEWKMNPTNCSTELRLLYQLVFLLEAHTCIPEN NGGAGCVCHLLMDDVVSADNYTLDLWAGQQLLWKGSFKPSEHVKPRAPGNL TVHTNVSDTLLLTWSNPYPPDNYLYNHLTYAVNIWSENDPADFRIYNVITYLEPS LRIAASLTKSGISYRARVRAWAQCYNNTTWSEWSPSTKWHNSYREPFEQH
Research Area	Immunology
Source	Mammalian cell
Target Names	IL4R
Expression Region	26-232aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	24.4 kDa
Protein Length	Extracellular Domain

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Endotoxin

Less than 1.0 EU/ μ g as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.