




# Recombinant Human Interleukin-2 (IL2) (Active)

<b>Product Code</b>	CSB-AP004381HU
<b>Abbreviation</b>	Recombinant Human IL2 protein (Active)
<b>Uniprot No.</b>	P60568
<b>Form</b>	Lyophilized powder
<b>Storage Buffer</b>	Lyophilized from a 0.2 μm Filtered 10 mM sodium citrate, pH 4.0
<b>Product Type</b>	Interleukin
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Biological Activity</b>	The specific activity of recombinant Human IL-2 as measured in a cell proliferation assay using CTLL-2 mouse cytotoxic T cells is $\geq 1 \times 10^7$ IU/mg.
<b>Purity</b>	Greater than 95% as determined by SDS-PAGE.
<b>Sequence</b>	APTSSSTKKTQLQLEHLLLDLQMILNGINNYKNPKLTRMLTFKGYMPKKATELK HLQCLEEELKPLEEVNLNAQSKNFHLRPRDLISNINVIVLELKGSETTFMCEYAD ETATIVEFLNRWITFCQSIISTLT
<b>Research Area</b>	Immunology
<b>Source</b>	E.coli
<b>Target Names</b>	IL2
<b>Expression Region</b>	21-153aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag-Free
<b>Mol. Weight</b>	15.5 kDa
<b>Protein Length</b>	Full Length of Mature Protein

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

<b>Endotoxin</b>	Less than 0.01 EU/μg as determined by LAL method.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.