




# Recombinant Human Interleukin-7 (IL7) (Active)

<b>Product Code</b>	CSB-AP004491HU
<b>Abbreviation</b>	Recombinant Human IL7 protein (Active)
<b>Uniprot No.</b>	P13232
<b>Form</b>	Lyophilized powder
<b>Storage Buffer</b>	Lyophilized from a 0.2 μm filtered 20 mM PB, 150 mM NaCl, pH 7.4
<b>Product Type</b>	Interleukin
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Biological Activity</b>	The ED50 as determined in a cell proliferation assay using PHA-activated human peripheral blood mononuclear cell (PBMC) is 50-300 pg/ml.
<b>Purity</b>	Greater than 95% as determined by SDS-PAGE.
<b>Sequence</b>	DCDIEGKDGKQYESVLMVSIQQLDSMKEIGSNCLNNEFNFFKRHICDANKEG MFLFRAARKLRQFLKMNSTGDFDLHLLKVSEGTILLNCTGQVKGRKPAALGE AQPTKSLEENKSLKEQKLLNDLCLFLKRLLEIKTCWNKILMGTKEH
<b>Research Area</b>	Immunology
<b>Source</b>	Mammalian cell
<b>Target Names</b>	IL7
<b>Expression Region</b>	26-177aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal 6xHis-tagged
<b>Mol. Weight</b>	18.4 kDa
<b>Protein Length</b>	Full Length of Mature Protein

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

<b>Endotoxin</b>	Less than 0.01 EU/μg as determined by LAL method.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.