




Recombinant Human Interleukin-31 (IL31) (Active)

Product Code	CSB-AP004631HU
Uniprot No.	Q6EBC2
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 μ m filtered 1xPBS, pH 7.4.
Product Type	Interleukin
Immunogen Species	Homo sapiens (Human)
Biological Activity	The ED50 as determined by inducing STAT3 activation in U87 MG human glioblastoma/astrocytoma cells is less than 45 ng/mL.
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	SHTLPVRLLRPSDDVQKIVEELQSLSKMLLKDVVEEEKGVLVSQNYTLPCLSPDAQPPNNIHSPAIRAYLKTIRQLDNKSVIDEIIEHLDKLIFQDAPETNISVPTDTHECKRFILTISQQFSECMDLALKSLTSGAQQATT
Research Area	Immunology
Source	E.coli
Target Names	IL31
Expression Region	24-164aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	15.8 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Endotoxin	Less than 0.01 EU/ μ g as determined by LAL method.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.