




Recombinant Mouse Interleukin-17A (Il17a), partial (Active)

Product Code	CSB-AP004781MO
Abbreviation	Recombinant Mouse Il17a protein, partial (Active)
Uniprot No.	Q62386
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered solution of PBS,PH7.4
Product Type	Interleukin
Immunogen Species	Mus musculus (Mouse)
Biological Activity	The ED50 as determined by its ability to induce IL-6 secretion by NIH?3T3 mouse embryonic fibroblast cells is 0.25-1.25 ng/ml.
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	TVKAAAIIPQSSACPNTAEKDFLQNVKVNLIKVFNSLGAKVSSRRPSDYLNIRSTS PWTILHRNEDPDRYPSVIWEAQCRHQRCVNAEGKLDHHMNSVLIQQEILVLKR EPESCPFTFRVEKMLVGVGCTCVASIVRQAA
Research Area	Immunology
Source	Mammalian cell
Target Names	Il17a
Expression Region	22-158aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	16.2 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Endotoxin

Less than 0.01 EU/µg as determined by LAL method.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.