




Recombinant Mouse Interleukin-4 (IL4) (Active)

Product Code	CSB-AP004831MO
Abbreviation	Recombinant Mouse IL4 protein (Active)
Uniprot No.	P07750
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 μ m filtered 20 mM PB, 150 mM NaCl, pH 7.4
Product Type	Interleukin
Immunogen Species	Mus musculus (Mouse)
Biological Activity	The ED50 as determined in a cell proliferation assay using M-NFS-60 mouse lymphoblast cells is 0.035 ng/ml.
Purity	\geq 95% as determined by SDS-PAGE.
Sequence	HIHGCDKNHLREIIGILNEVTGEGTPCTEMDVPNVLTKNTTESELVCRASKVL RIFYLKHGKTPCLKKNSSVLMELQRLFRAFRCCLDSSISCTMNESKSTSLKDFLE SLKSIMQMDYS
Research Area	Immunology
Source	Mammalian cell
Target Names	IL4
Expression Region	21-140aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	14.6 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Endotoxin	Less than 1.0 EU/ μ g as determined by LAL method.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.