




Recombinant Human Tumor necrosis factor ligand superfamily member 4 (TNFSF4), partial (Active)

Product Code	CSB-AP005261HU
Uniprot No.	P23510
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 μ m filtered 20 mM PB, 150 mM NaCl, pH 7.4
Product Type	Immune Checkpoint
Immunogen Species	Homo sapiens (Human)
Biological Activity	The ED50 as determined by its ability to bind Human OX40 in functional ELISA is less than 20 ug/ml.
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	QVSHRYPRIQSIKVFTEYKKEKGFILTSQKEDEIMKVQNNNSVIINCDGFYLISLK GYFSQEVNISLHYQKDEEPLFQLKKVRSVNSLMVASLTYKDKVYLNVTDDNTS LDDFHVNGGELILIHQNPGEFCVL
Research Area	Cancer
Source	Mammalian cell
Target Names	TNFSF4
Expression Region	51-183aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	16.3 kDa
Protein Length	Extracellular Domain

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Endotoxin	Less than 1.0 EU/ μ g as determined by LAL method.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.