




# Recombinant Human Hepatitis A virus cellular receptor 2 (HAVCR2), partial (Active)

<b>Product Code</b>	CSB-AP005281HU
<b>Uniprot No.</b>	AAL65157.1
<b>Storage Buffer</b>	Lyophilized from a 0.2 μm filtered 1xPBS, pH 7.4
<b>Product Type</b>	Immune Checkpoint
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sensitivity</b>	The ED50 as determined by its ability to bind Human Galectin 9 in functional ELISA is less than 15 ug/ml.
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	SEVEYRAEVLGQNAAYLPCFYTPAAPGNLVPVCWVGKACPVFECGNVVLRTDE RDVNYWTSRYWLNQDFRKGDVSLTIENVTLADSGIYCCRIQIPGIMNDEKFNK LVIKPAKVTPAPTLQRDFTAAFPRLTTRGHGPAETQTLGSLPDINLTQISTLAN ELRDSRLANDLRDSGATIR
<b>Research Area</b>	Immunology
<b>Source</b>	Mammalian cell
<b>Target Names</b>	HAVCR2
<b>Expression Region</b>	22-200aa
<b>Tag Info</b>	C-terminal FC-Avi-tagged
<b>Mol. Weight</b>	48.8 kDa
<b>Protein Length</b>	Partial

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

<b>Typical</b>	Less than 1.0 EU/μg as determined by LAL method.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.