




Recombinant Human Carbonic anhydrase 7 (CA7) (Active)

Product Code	CSB-AP005411HU
Uniprot No.	P43166
Form	Liquid
Storage Buffer	0.2 μ m filtered 20 mM Tris-HCl, 150 mM NaCl, pH 8.0
Product Type	Other
Immunogen Species	Homo sapiens (Human)
Biological Activity	The esterase activity is determined to be greater than 1000 pmol/min/ug
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	MTGHHGWGYGQDDGPSHWHKLYPIAQGDRQSPINIISSQAVYSPSLQPLELS YEACMSLSITNNGHSVQVDFNDSDDRTVVTGGPLEGPYRLKQFHFHWGKKH DVGSEHTVDGKSFPSELHLVHWNACKYSTFGEAASAPDGLAVVGVFLETGDE HPSMNRLTDALYMVRFKGTAKQFSCFNPKCLLPASRHYWTYPGSLTTPPLSE SVTWIVLREPICISERQMGKFRSLLFTSEDDERIHMVNNFRPPQPLKGRVVKAS FRA
Research Area	Cell Biology
Source	E.coli
Target Names	CA7
Expression Region	1-264aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	30.72 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Endotoxin**

Less than 1.0 EU/ μ g as determined by LAL method.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.