




# Recombinant Human Ephrin type-B receptor 1 (EPHB1), partial (Active)

<b>Product Code</b>	CSB-AP005501HU
<b>Uniprot No.</b>	P54762
<b>Form</b>	Lyophilized powder
<b>Storage Buffer</b>	Lyophilized from a 0.2 $\mu$ m Filtered 20 mM Tris-HCl, 150 mM NaCl, pH 8.0
<b>Product Type</b>	Other
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Biological Activity</b>	The ED50 as determined by its ability to binding Human EFNB2 used funtional ELISA is less than 100 ug/ml.
<b>Purity</b>	$\geq$ 95% as determined by SDS-PAGE.
<b>Sequence</b>	SRKRAYSKEAVYSDKLQHYSTGRGSPGMKIYIDPFTYEDPNEAVREFAKEIDV SFVKIEEVIGAGEFGEVYKGRLLKLPKREIYVAIKTLKAGYSEKQRRDFLSEASI MGQFDHPNIIRLEGVVTKSRPVMIIITEFMENGALDSFLRQNDGQFTVIQLVGML RGIAAGMKYLAEMNYVHRDLAARNILVNSNLVCKVSDFGLSRYLQDDTSDPTY TSSLGGKIPVRWTAPEAIAYRKFTSASDVWSYGIVMWEVMSFGERPYWDMSN QDVINAIEQDYRLPPPMDCPAALHQLMLDCWQKDRNSRPRFAEIVNTLDKMIR NPASLKT VATITAVPSQPLLD RSIPDFTAFTTVDDWLSAIK MVQYRDSFLTAGFT SLQLVTQMTSEDLLRIGITLAGHQKKILNSIHS MRVQISQSPTAMA
<b>Research Area</b>	Signal Transduction
<b>Source</b>	Mammalian cell
<b>Target Names</b>	EPHB1
<b>Expression Region</b>	564-984aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal 6xHis-tagged
<b>Mol. Weight</b>	48.8 kDa
<b>Protein Length</b>	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Endotoxin**

Less than 1.0 EU/μg as determined by LAL method.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.