




Recombinant Human Parathyroid hormone (PTH) (Active)

Product Code	CSB-AP005621HU
Uniprot No.	P01270
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 μ m Filtered 10 mM HAc-NaAc, 150 mM NaCl, 5% Mannitol, pH 4.0
Product Type	Other
Immunogen Species	Homo sapiens (Human)
Biological Activity	The ED50 as determined by its ability to induce cAMP accumulation in MC3T3?E1 mouse preosteoblast cells is less than 0.1 ug/ml
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	SVSEIQLMHNLGKHLNSMERVEWLRKKLQDVHNFVALGAPLAPRDAGSQRPR KKEDNVLVESHEKSLGEADKADVNVLTAKASQ
Research Area	Signal Transduction
Source	E.coli
Target Names	PTH
Expression Region	32-115aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	9 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Endotoxin	Less than 1.0 EU/ μ g as determined by LAL method.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.