



ALOX15 Antibody, Biotin conjugated

Product Code	CSB-PA001621LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16050
Immunogen	Recombinant Human Arachidonate 15-lipoxygenase protein (1-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arachidonate 15-lipoxygenase (15-LOX) (15-LOX-1) (EC 1.13.11.33) (12/15-lipoxygenase) (Arachidonate 12-lipoxygenase, leukocyte-type) (12-LOX) (EC 1.13.11.31) (Arachidonate omega-6 lipoxygenase), ALOX15, LOG15
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ALOX15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.