



APC2 Antibody, Biotin conjugated

Product Code	CSB-PA001895LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95996
Immunogen	Recombinant Human Adenomatous polyposis coli protein 2 protein (2011-2221AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenomatous polyposis coli protein 2 (Adenomatous polyposis coli protein-like) (APC-like), APC2, APCL
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	APC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.