



ERAP2 Antibody, Biotin conjugated

Product Code	CSB-PA007761LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P179
Immunogen	Recombinant Human Endoplasmic reticulum aminopeptidase 2 protein (406-522AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoplasmic reticulum aminopeptidase 2 (EC 3.4.11.-) (Leukocyte-derived arginine aminopeptidase) (L-RAP), ERAP2, LRAP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ERAP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.