



ETV2 Antibody, Biotin conjugated

Product Code	CSB-PA007855LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00321
Immunogen	Recombinant Human ETS translocation variant 2 protein (74-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ETS translocation variant 2 (Ets-related protein 71), ETV2, ER71 ETSRP71
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ETV2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.